



## Feed applications for trinamiX PAL One

Know exactly what is in your feed materials –  
as ingredients or finished feed

# trinamiX

A brand of  
BASF – We create chemistry

Model Parameter	Finished Feed			Feed Ingredients					Forages			
	 Ruminant	 Swine	 Poultry	 Cereals	 Oilseeds & Expeller Meals	 Extraction Meals	 By-products	 Animal Proteins	 Grain Silage	 Grass Silage	 Hay	 Fresh Grass
Moisture	✓	✓	✓	✓	✓	✓	✓	✓				
Dry Matter									✓	✓	✓	✓
Crude Protein	✓	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓
Protein									✓			
Fat (Ether Extract)	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Starch	✓	✓	✓	✓	✓	✓	✓		✓			
Sugar	✓	✓	✓	✓	✓	✓	✓					
Crude Fibre	✓	✓	✓	✓	✓	✓	✓	✓				
Fat (Acid Hydrolysis)	✓	✓	✓	✓	✓	✓	✓	✓				
NDF	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓
ADF	✓	✓	✓	✓	✓	✓	✓		✓	✓	✓	✓
Ash	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Amino Acids				✓	✓	✓	✓	✓				
Digestible Energy	✓	✓							✓	✓	✓	✓
Metabolizable Energy	✓	✓							✓	✓	✓	✓
Net Energy Lactation	✓								✓	✓	✓	✓
Gross Energy		✓										
Net Energy		✓										
AMEn			✓	✓								
D-Value									✓	✓	✓	✓
WSC										✓	✓	✓